

钢制切割模具

型号：C0024

IDM[®]
instruments

钢制切割模具主要用于切割塑料，纸张，纺织品和橡胶等试样。它可以与 IDM 的切割机一起使用。

产品应用：

- 纸板
- 海绵
- 纸张
- 塑料
- 橡胶

产品特点：

- 可以反复磨快
- 比其他材质的模具使用寿命长

参考标准：

- BS 2782-3
- AS 1301.404
- AS 1683.11
- AS 1180.2
- AS 2282.7
- ASTM D412
- ASTM D470
- ASTM D624
- ASTM D1004
- ASTM 1660.2
- ASTM D1822
- ASTM D3574
- ASTM D638
- AS 1145
- ISO 34-1
- ISO 37
- ISO 527
- ISO 1798
- ISO 6259
- K6251
- K6301
- AS/NZS 1660.2.3
- ASTM D882

选配辅具：

- 定制模具-客户要求
- 切割模具适配器
- 切割机
- 证书

模具尺寸：

- 参考以下图纸（单位为 mm）



C0022

C0050

S0003




澳大利亚 IDM 仪器公司.北京代表处
北京蓝瑞泰克科技有限公司
电话：010-67086989
传真：010-67085997
邮箱：idm@laryee.com
网站：www.idminstruments.cn

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 – Type A

- ISO 37-Type 2
- ISO 527-2-Type 5A
- AS 1145-Type 5A
- K6251-Type 6
- AS/NZS 1660.2.3
- DIN53504




C0024 – Type B

- ASTM D412-C
- ISO 6259-3-Type 2
- K6251-Type 5
- ISO 37 – Type 1
- ISO 527-3 – Type 5
- BS 2782-3 (Method 320A)
- AS 1145.3 – Type 5
- DIN53504
- ASTM D882



C0024 – Type C

- ASTM D412-D



C0024 – Type D


- JIS K6301-Type 3
- K6251- Type 3
- ISO 37 – Type 1A

澳大利亚 IDM 仪器公司.北京代表处
北京蓝瑞泰克科技有限公司
电话: 010-67086989
传真: 010-67085997
邮箱: idm@laryee.com
网站: www.idminstruments.cn

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 - Type E

- ASTM D638-02a-Type 4
- ASTM D6893




C0024 - Type F

- JIS K6301-Type 1



C0024 - Type G

- JIS K6301-Type 2




C0024 - Type H

- JIS K6301-Type 4

钢制切割刀具


型号: C0024

IDM®
instruments




C0024 – Type I

- ISO 37 – Type 3
- AS/NZS 1660.2.3
- DIN 53504




C0024 – Type J

- ISO 37 – Type 4
- ISO 527-2 – Type 5B



C0024 – Type K

- K6301-B
- ASTM D624 – Type C
- ASTM D1004-03
- DIN ISO34-1




C0024 – Type L

- ASTM D638-02a-Type 1

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 – Type M

- ASTM D638-02a-Type 2




C0024 – Type N

- ASTM D638-02a-Type 3



C0024 – Type O

- ASTM D638-02a-Type 5




C0024 – Type P

- DIN ISO 34-1
- K6301-Type A
- ASTM D624 – Type B

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 – Type Q

- ASTM D1822-Type L




C0024 – Type R

- ASTM D1822-Type S



C0024 – Type S

- ASTM D3574 – Type E




C0024 – Type T

- ASTM D3574 – Type F
- DIN ISO34-1

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 - Type U

- ISO 527-2 - Type 1A




C0024 - Type V

- ISO 527-2 - Type 1B
- ISO 527-3 - Type 1B
- BS 2782-3 (Method 320B)
- BS 2782-10 - Type 1
- ISO 6259-3 Type 1




C0024 - Type W

- ISO 527-2 - Type 1BA



C0024 - Type X

- ISO 527-2 - Type 1BB



C0024 - Type Y


- ISO 527-3 - Type 2
- AS 1145.2 - Type 2

澳大利亚 IDM 仪器公司.北京代表处
北京蓝瑞泰克科技有限公司
电话: 010-67086989
传真: 010-67085997
邮箱: idm@laryee.com
网站: www.idminstruments.cn

钢制切割刀具


型号: C0024

IDM[®]
instruments




C0024 – Type Z

- ISO 527-3 – Type 4




C0024 – Type Aa

- ASTM D412 – Type A



C0024 – Type Ab

- ASTM D412 – Type B




C0024 – Type Ac

- ASTM D412 – Type E


钢制切割刀具

型号: C0024



C0024 – Type Ad

• ASTM D412 – Type F



C0024 – Type Ae

• ISO 1798