

Dry Bath Aging Block

Model: D0001

a measurable difference...

IDM[®]
instruments

To be used for determining the long-term effects of elevated temperature on various materials like rubber, plastics, textile etc. Samples are placed in glass tubes. The 24 holes for insertion of the tubes are 38mm by 280mm.

Applications:

- Rubber
- Plastic
- Textiles

Features:

- Set of glassware
- Temperature controller: Ambient – 300 °C

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- N/A

Options:

- N/A

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ
(please specify when ordering)

Dimensions:

- **H:** 500mm
- **W (OD):** 565mm
- **Control Box:** 150mm
- **Weight:** 90kg

