

Polarising Box

Model: P0008

a measurable difference...

IDM[®]
instruments

The Polarising Box is used to view the stress points created during manufacture of clear plastic and glass products.

The samples can also be extruded film.

The Polarising Box has two open sides so long items can be passed through.

Applications:

- Plastic films
- Glass

Features:

- Polarising Filter mounted at the top for viewing
- Light Source
- Fluorescent Tubes, mounted at base with a diffuser as the sample support panel

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- N/A

Options:

- Larger Polarising Box (Consult Factory)

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ

(please specify when ordering)

Dimensions:

- **H:** 400mm
- **W:** 400mm
- **D:** 400mm
- **Weight:** 7kg

