

Colour Fastness Tester

Model: C0010

a measurable difference...

IDM[®]
instruments

This apparatus is one of a series for determining the colour fastness of textiles to various agencies. It determines the colourfastness of textiles to light using an artificial light source.

Applications:

- Fabrics
- Fibres
- Yarn
- NOT polychromatic coloured textiles

Features:

- 500 Watt lamp
- 12 Sample holders
- Black panel thermometer
- Hour Timer
- Countdown Timer

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- AS 2001.4.21

Options:

- Centering Device

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ or
110 VAC @ 60 HZ

(please specify when ordering)

Dimensions:

- **H:** 400mm
- **W:** 650mm
- **D:** 500mm
- **Weight:** 10kg

