

Ningbo Abrasion Type Tester

Model: M0013

a measurable difference...

IDM[®]
instruments

A Martindale Abrasion & Pilling Tester used for testing the wearability and pilling properties of wool, wool blended, woven, knitted fabrics.

Application:

- Textiles

Features:

- Two trajectories and a reciprocating linear track.
- Capable of round tracking movement.
- Touch screen operation
- LCD Display
- Six operational preset procedures (ie. Repeat single procedure or three velocities can be set)
- Automatic stop

Specifications:

- 8 Stations
- Counter Display:
 - Preset counter display 0-999999
 - Cumulative counter display 0-999999
- Testing Pressure:
 - A: Cloth specimen clasper 155g ± 1g
 - B: Cloth specimen sample drag weight 260g ± 1g
- Loading drag weight 2385g ± 10g
- Abrading Table Diameter 120

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- ASTM D4966/4970
- BS3424/5690
- BS EN 388/530
- DIN53863/53865
- GB/T4802.2
- ISO5470
- IWS TM112/196
- JIS L1096

Connections:

- **Electrical:** 220/240 VAC @ 50 HZ



Dimensions:

- **H:** 365 mm
- **W:** 820 mm
- **D:** 605 mm
- **Weight:** 168kg