

Digital Air Permeability Tester

Models: A0004

a measurable difference...

IDM[®]
instruments

The Digital Air Permeability Tester determines the air leakage of plastic films. It is manufactured to enable samples to be inserted at the front, with a sample opening of 5mm as well as a 45° lead in on sample opening. It features two safety activation switches for clamping, so that operators are using both hands to ensure optimum safety.

Applications:

- Plastic

Features:

- Sample opening: 5mm
- Low pressure air supply gauge
- Digital display
- 2 x safety activation buttons
- Test area: 250mm x 120mm
- Air flow meter - flow range: 0-50 LPM for air @ 1 atmospheric, 20°C (sensed before sample)

Connections:

- **Air supply:** 60 psi
- **Electrical:** 220/240 VAC @ 50 HZ or 110 VAC @ 60 HZ
(please specify when ordering)

Dimensions:

- **H:** 400mm
- **W:** 330mm
- **D:** 360mm
- **Weight:** 40kg

