

Cutting Dies – Ruler Type

Model: C0025

a measurable difference...

IDM[®]

instruments

The Ruler Type Cutting Dies are used for cutting of plastic film, paper and rubber samples used for tensile and tear testing (dumbbell, trouser, etc). They are manufactured with a ply backed knife-edge, and can be manufactured to suit various international test standards.

Applications:

- Plastic Film
- Paper
- Rubber

Features:

- Low volume work

Benefits:

- Easy to use
- Fast results
- Accurate

Standards:

- BS2783
- AS1683.11
- AS2282.7
- ASTM D470
- ASTM D638
- ASTM 1660.2
- ASTM D3574
- AS1301.404
- AS1180.2
- ASTM D412
- ASTM D624
- ASTM D1004
- ASTM D1822

Options:

- Cutting Press

Connections:

- N/A

Dimensions:

- As per requirements



Can be used with the following Cutting Presses

