

# 刀模-尺型

型号: C0025

**IDM**<sup>®</sup>  
instruments

此模具用于切割塑料薄膜、纸张、橡胶样品（哑铃型等）做拉伸、撕裂测试。可用刀来进行手工切割，也可配合各款切压机使用。

## 产品应用:

- 塑料薄膜
- 纸张
- 橡胶

## 产品特点:

- 性价比高

## 产品优势:

- 使用简单
- 测试快速
- 结果精确

## 参考标准:

- BS 2783
- AS 1683. 11
- AS 2282. 7
- ASTM D470
- ASTM D638
- ASTM 1660. 2
- ASTM D3574
- AS 1301. 404
- AS1180. 2
- ASTM D412
- ASTM D624
- ASTM D1004
- ASTM D1822

## 选配附件:

- 切割器

