Cutting Press-3 Tonne Arbour a measurable difference... Model: C0022 - 3

The handle operated manual Arbour Press is used in conjunction with cutting dies to produce testing samples. An optional cutting die adaptor can be used when using the Arbour Press with the Tool Steel cutting dies for fast and repetitious sample preparation.

Applications:

- Paper
- Plastic
- Rubber
- Textiles

Features:

- Ease of operation
- Rack and Pinion
- Table size: 250 mm x 300 mm
- Pressure Foot: 125mm Ø

Benefits:

- · Easy to use
- Fast results

Standards:

• N/A

Options:

- Cutting Die Adaptor
- Cutting Dies

Connections:

• N/A

Dimensions:

- **H:** 650mm
- W: 350mm
- **D:** 450mm
- Weight: 73kg



Cutting Die Adaptor



Appearance & specifications listed are subject to change without notice. Copyright © 2011 IDM Instruments Pty Ltd. All Rights Reserved.







Ph: +61 3 9708 6885 Fax: +61 3 9708 6770 Email: idm@idminstruments.com.au Web: www.idminstruments.com.au