

Crockmeter

Model: C0019

a measurable difference...

IDM[®]
instruments

The specimen is clamped to the base plate, where a peg is moved back and forth over the test specimen to determine colour fastness from dry or wet tests.

Applications:

- Textiles
- Carpets



Features:

- Testing Times: 1-9999 (LED)
- Rubbing Travel: 104mm ± 3mm
- Rubbing Pressure: 9N
- Rubbing Speed: 60r/min
- Rubbing Head: Round 16mm Ø
Rectangle 19mm x 25mm

Standards:

- AS 2001.4.3
- AS / NZ 2111.19.1
- AATCC - Test method 8-1972
- BS 1006
- GB/T3920
- ISO 105
- IWS TM165
- JIS K6328

Benefits:

- Easy to use
- Fast results
- Accurate

Connections:

- AC 220V

Dimensions:

- **H:** 280mm
- **L:** 770mm
- **D:** 250mm
- **Weight:** 25kg